

## EthosButton CL

UHMWPE Nonabsorbable Suture Loop

---

HNM Medical has developed the **EthosButton CL** a product that is projected to have comparable strength to similar devices indicated for anterior cruciate ligament (ACL) femoral cortical fixation repair. The **EthosButton CL** consists of an ultra-high molecular weight polyethylene (UHMWPE) loop, combined with a titanium plate (button).

### The EthosButton CL Advantage

We analyzed comparative loop and button strength using cyclical tension and static (ultimate) load testing. Cyclical tension testing was performed using an Electromechanical Dynamic Testing System at a 1000N force. Ultimate load failure was determined using an Electromechanical Dynamic Testing System.

Available Log Sizes: 15mm, 20mm, 25mm, 30mm, 35mm





For more information please contact us at [info@hnmmedical.com](mailto:info@hnmmedical.com) or 1.866.291.8489  
to view our complete portafolio visit us at: [www.hnmmedical.com](http://www.hnmmedical.com)