



**NAOS First
MPJ Hemi
Implant System**

PRODUCT INFORMATION





Contents

| | |
|----------------------|---|
| Introduction | 3 |
| Features | 3 |
| Instrumentation | 4 |
| Ordering Information | 5 |



HNMTOTAL RECON is committed in providing surgeons and health care professionals with advanced surgical solutions for effective and time-saving treatment of their patients. Our most important task is meeting surgeons' needs. We accomplish this by creating key partnerships and driving innovative product development.

📍 20855 Northeast 16th Avenue, Suite C15
Miami, Florida 33179

☎ 1.866.291.8498 ✉ info@hnmtotal.com
📞 1.888.608.9788 🌐 www.hnmtotal.com


NAOS FIRST MPJ HEMI

IMPLANT SYSTEM

The NAOS First MPJ Hemi Implant is simple to implant, removes the minimum amount of bone possible, and preserves the anatomy of the first toe. The NAOS implant is an excellent alternative to fusion, offering patients pain relief and restored range of motion.

The implant offers many benefits over newer, unproven implants designed for the metatarsal head. Multiple clinical studies have shown the implant to be successful with no cartilage on the metatarsal head, countering the belief that you must treat the side of the joint with visible cartilage damage. The base of the phalanx is subject to much lower forces than the metatarsal head (which can be subject to six times body weight when weight bearing). Additionally, the NAOS First MPJ Hemi Implant surgical procedure can be combined with distal metatarsal osteotomies for correction of accompanying pathologies.

FEATURES






INSTRUMENTATION

NAOS FIRST MPJ HEMI BASIC INSTRUMENT SET

HTR-01-730768



CONTENT

| Item # | Description | Pcs. |
|---------------|---|------|
| HTR-01-770186 | NAOS Punch | 1 |
| HTR-01-730209 | NAOS Template 5-Prong | 1 |
| HTR-01-220335 | NAOS Trial SM | 1 |
| HTR-01-710480 | NAOS Trial MS | 1 |
| HTR-01-220536 | NAOS Trial MD | 1 |
| HTR-01-490659 | NAOS Trial ML | 1 |
| HTR-01-220737 | NAOS Trial LG | 1 |
| HTR-01-510812 | NAOS Chisel | 1 |
| HTR-01-520359 | Impactor | 1 |
| HTR-01-520256 | Impactor Tip SM | 1 |
| HTR-01-520957 | Impactor Tip LG | 1 |
| HTR-01-821000 | Clamp | 1 |
| HTR-01-800160 | Tray for NAOS First MPJ Hemi Basic Instrument Set | 1 |

ORDERING INFORMATION

NON-POROUS COATED



| Toe Joint | Size | Cobalt Chrome |
|-----------|--------|---------------|
| S | 17mm | HTR-01-000760 |
| MS | 18.5mm | HTR-01-700833 |
| M | 20mm | HTR-01-000961 |
| ML | 21.5mm | HTR-01-491058 |
| L | 23mm | HTR-01-000162 |

AVAILABLE IN TITANIUM
FOR PATIENTS WITH METAL
SENSITIVITY.

POROUS COATED



| Toe Joint | Size | Cobalt Chrome |
|-----------|--------|---------------|
| S | 17mm | HTR-01-040212 |
| MS | 18.5mm | HTR-01-700334 |
| M | 20mm | HTR-01-040413 |
| ML | 21.5mm | HTR-01-490560 |
| L | 23mm | HTR-01-040614 |

HNMTOTAL

TOTAL RECON

📍 20855 NE 16 Ave Suite C15, Miami, FL 33179

☎ 1.866.291.8498

✉ info@hnmtotal.com

📠 1.888.608.9788

🌐 www.hnmtotal.com