




HNM

EMPOWERING INNOVATION

MULTIFUNCTIONAL TABLE MOUNTED RETRACTOR SYSTEM

MULTIFUNCTIONAL TABLE MOUNTED RETRACTOR SYSTEM

All parts are self-locking. This reliable retractor system provides the surgeon with maximum surgical exposure.






Oval Frame

Small	HNM-4-600-011	203x311mm
Medium	HNM-4-600-012	267x368mm
Large	HNM-4-600-013	267x470mm



Round Frame

Small	HNM-4-600-021	216mm D
Medium	HNM-4-600-022	267mm D
Large	HNM-4-600-023	317mm D



Oval Segmented Frame

Small	HNM-4-600-031	216x343mm
Medium	HNM-4-600-032	267x495mm



Round Segmented Frame

Small	HNM-4-600-041	216mm D
Medium	HNM-4-600-042	267mm D



Vaginal Oval Frame

HNM-4-600-051	203x266mm
---------------	-----------



Segments (pair with 2 pieces)

HNM-4-600-061	102mm
HNM-4-600-062	152mm
HNM-4-600-063	203mm


Balfour retractor blade, 220mm length

HNM-4-600-201	68x85mm
HNM-4-600-202	70x75mm




Kelly retractor blade, 180mm length

HNM-4-600-211	63x51mm
HNM-4-600-212	51x51mm
HNM-4-600-213	75x65mm
HNM-4-600-214	76x51mm
HNM-4-600-215	102x51mm
HNM-4-600-216	127x51mm
HNM-4-600-217	154x51mm
HNM-4-600-218	190x50mm




Harrington heartshaped retractor blade, 200mm length

HNM-4-600-221	145x64mm
---------------	----------




Rectal retractor blade, 167mm length

HNM-4-600-311	75mm
HNM-4-600-312	95mm



Hook retractor blade, 190mm length

HNM-4-600-251	25.4X25.4mm
HNM-4-600-252	51X25.4mm
HNM-4-600-253	102X25.4mm




Malleable retractor blade, 302mm length

HNM-4-600-231	150X38mm
HNM-4-600-232	150X51mm
HNM-4-600-233	150X76mm



Deaver retractor blade, 215mm length

HNM-4-600-261	89x25mm
HNM-4-600-262	127x25mm



Gelpi hook

HNM-4-600-241	165mm
---------------	-------






Vaginal retractor blade #1,
180mm length

HNM-4-600-271 83x25mm




Vaginal retractor blade #2,
215mm length

HNM-4-600-272 114x25mm




Kocher retractor blade, 180mm length

HNM-4-600-281	38x51mm
HNM-4-600-282	40x25mm
HNM-4-600-283	65x25mm
HNM-4-600-284	122x20mm
HNM-4-600-285	122x25mm




Richardson retractor blade, 180mm length

HNM-4-600-291	25X19mm
HNM-4-600-292	51X19mm
HNM-4-600-293	34X25mm



Meyerdig serrated retractor blade,
with cross bar, 175mm length

HNM-4-600-321	75X26mm
HNM-4-600-322	78X25mm




Meyerdig serrated retractor blade,
180mm length

HNM-4-600-331	62X17mm
HNM-4-600-332	81X25mm
HNM-4-600-333	92X50mm



Skin rake retractor, 188mm length

HNM-4-600-343	3 teeth 17mm wide
HNM-4-600-344	4 teeth 19mm wide
HNM-4-600-345	6 teeth 25mm wide



Kelly vein retractor blade, 180mm length

HNM-4-600-351	76x21mm
HNM-4-600-352	102x21mm
HNM-4-600-353	127x21mm
HNM-4-600-354	152x21mm
HNM-4-600-355	178x21mm
HNM-4-600-356	203x21mm



Table post clamping device
HNM-4-600-100



Horizontal bar
HNM-4-600-102



Horizontal bar with flexible head
HNM-4-600-103




Extension bar
HNM-4-600-101



Post coupler
HNM-4-600-104



Tilt ratchet mechanism
HNM-4-600-110



Standard ratchet
HNM-4-600-111



Fibreoptic light D3.5mm, 180mm length
HNM-4-798-060



TOLL FREE 1.866.291.8498
FAX 1.888.608.9788
WWW.HNMMEDICAL.COM
E-MAIL: INFO@HNMMEDICAL.COM
