



# Distraction Pin

Distraction Pin is designed for Anterior Cervical Fusion Procedures. Use with vertebral body distraction.

The device is available in varying tip lengths.

Each pin is manufactured from 416 SS, COND, T per ASTM A582

Part#	Tip Length (mm)	Description
HNM 909 - 12	12	Distraction pin
HNM 909 - 14	14	Distraction pin
HNM SS - 0088S	12	Distraction pin, sharp
HNM SS -0089S	14	Distraction pin, sharp



Customer Service and Ordering  
HNM Medical – USA  
20855 NE 16 Avenue Suite C15  
Miami, FL 33179  
T: 1.866.291.8498  
F: 1.888.608.9788  
E: [info@hnmmedical.com](mailto:info@hnmmedical.com)

