



HNM **SUPPLIES**

PRODUCT INFORMATION



GAUZE
COTTON

BANDAGE
TAPE

NON-WOVEN
DRAPE & BACK

UNDERPAD DIAPER
INFUSION SET

DISPOSABLE PRODUCT

CONTENTS

GAUZE

Sterile lap sponge	3
Non sterile lap sponge	4

COTTON

Cotton roll	5
-------------	---

BANDAGE

POP Bandage	
- 5cmx2.7m	6
- 10cmx2.7m	6
- 20cmx2.7m	6

TAPE

NON-Woven Tape	7
----------------	---

NON-WOVEN

Cap Doctor Cape	8
Coverall	9
AAMI Surgical gown	10
Transverse incision	11

DRAPE & BACK

General Pack	12
--------------	----

UNDERPAD DIAPER

Baby diaper	13
Adult diaper	13

INFUSION SET

Infusion set	14
--------------	----

DISPOSABLE PRODUCT

IV Canula	16
Needle for anesthesia	16
Feeding tube	17
Nebulizer Mask	17
Breathing circuit	18
Blood Bag	18
Adult urine bag	19
Autoclave Steril	19

GAUZE

STERILE LAP SPONGE

- 100% Cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 22, 20, 17, 15, 13, 12, 11 threads etc.
- With or without x-ray detectable/ x-ray film, with or without blue loop
- Pre-washed or non-washed
- 4ply, 6ply, 8ply, 12ply, 16ply, 24ply etc
- 20x20cm, 22.5x22.5cm, 30x30cm, 40x40cm, 45x45cm, 50x50cm etc
- 5's 2's, 1's in sterile pouch or double package
- Gamma, EO, steam

ITEM #	SPECS	PACKING	CARTON SIZE
	45cm x 45cm - 12ply	200	57 x 30 x 42cm
	40cm x 40cm- 12ply	200	48 x 30 x 38cm
	30cm x 30cm - 12ply	500	52 x 36 x 40cm
	30cm x 30cm - 12ply	100 pk	57 x 25 x 46cm
	45cm x 45cm - 8ply	400	57 x 42 x 42cm
	45cm x 45cm - 8ply	80 pk	60 x 28 x 50cm
	40cm x 40cm- 12ply	400	48 x 42 x 36cm
	30cm x 30cm - 8ply	600	57 x 36 x 45cm
	45cm x 45cm - 8ply	800	57 x 42 x 42cm
	45cm x 45cm - 8ply	200 pk	54 x 39 x 52cm
	40cm x 40cm- 12ply	800	48 x 42 x 38cm
	30cm x 30cm - 8ply	1200	57 x 40x 45cm
	30cm x 30cm - 8ply	200 pk	57 x 38 x 40cm



NON - STERILE LAP SPONGE

- 100% Cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 22, 20, 17, 15, 13, 12, 11 threads etc.
- With or without x-ray detectable/ x-ray film, with or without blue loop
- Pre-washed or non-washed
- 4ply, 6ply, 8ply, 12ply, 16ply, 24ply etc
- 20x20cm, 22.5x22.5cm, 30x30cm, 40x40cm, 45x45cm, 50x50cm etc
- 10's 25's, 50's 100's per pack (Non sterile)

ITEM #	SPECS	PACKING	CARTON SIZE
	45cm x 45cm - 12ply	200	57 x 30 x 42cm
	40cm x 40cm- 12ply	200	48 x 30 x 38cm
	30cm x 30cm - 12ply	500	52 x 36 x 40cm
	30cm x 30cm - 12ply	100 pk	57 x 25 x 46cm
	45cm x 45cm - 8ply	400	57 x 42 x 42cm
	45cm x 45cm - 8ply	80 pk	60 x 28 x 50cm
	40cm x 40cm- 12ply	400	48 x 42 x 36cm
	30cm x 30cm - 8ply	600	57 x 36 x 45cm
	45cm x 45cm - 8ply	800	57 x 42 x 42cm
	45cm x 45cm - 8ply	200 pk	54 x 39 x 52cm
	40cm x 40cm- 12ply	800	48 x 42 x 38cm
	30cm x 30cm - 8ply	1200	57 x 40x 45cm
	30cm x 30cm - 8ply	200 pk	57 x 38 x 40cm



COTTON ROLL

- Material: 100% pure cotton
- Packing: 1roll/blue kraft paper or polybag
- It is suite for medical and daily use
- Type: Ordinary precut

ITEM #	SPECS	PACKING	CARTON SIZE
	25 gr/roll	500 rolls/ctn	56 x 36 x 56cm
	40 gr/roll	400 rolls/ctn	56 x 37 x 56cm
	50 gr/roll	300 rolls/ctn	61 x 37 x 61cm
	80 gr/roll	200 rolls/ctn	61 x 31 x 61cm
	100 gr/roll	200 rolls/ctn	61 x 31 x 61cm
	125 gr/roll	100 rolls/ctn	61 x 36 x 36cm
	200 gr/roll	50 rolls/ctn	41 x 41 x 41cm
	250 gr/roll	50 rolls/ctn	41 x 41 x 41cm
	400gr/roll	40 rolls/ ctn	55 x 31 x 36cm
	454 gr/roll	40 rolls/ctn	61 x 37 x 46cm
	500 gr/roll	20 rolls/ctn	61 x 38 x 48cm
	1000 gr/roll	20 rolls/ctn	66 x 34 x 52cm



POP BANDAGE

- Material: cotton + plaster of paris
- Setting time: 2-3 minutes
- Width: 5cm, 7.5cm, 10cm, 12.5cm, 15cm, 30cm, etc.
- Length: 10m, 10 yards, 5m, 5 yards, 4m, 3m, 2.7m etc.
- Packing: 1roll/ individually packed

ITEM #	SIZE	PACKING	CARTON SIZE
	5cm x 2.7m	240 rolls/ctn	57 x 33 x 26cm
	7.5cm x 2.7m	240 rolls/ctn	57 x 33 x 26cm
	10cm x 2.7m	120 rolls/ctn	57 x 33 x 26cm
	12.5cm x 2.7m	120 rolls/ctn	57 x 33 x 26cm
	15cm x 2.7m	120 rolls/ctn	57 x 33 x 26cm
	20cm x 2.7m	60 rolls/ctn	57 x 33 x 26cm



NON - WOVEN TAPE

- Material: Non-woven
- Color: Skin o white color
- Glue: hot melt or acrilyc acid glue
- Simple packing / plastic tin / with dispenser
- Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
- Length: 10m, 10 yards, 5m, 5yards, etc.

ITEM #	SIZE	CARTON SIZE	PACKING
	1.25cm x 5yds	24 x 23.5 x 28.5	57 x 33 x 26cm
	2.5cm x 5 yds	24 x 23.5 x 28.5	57 x 33 x 26cm
	5cm x 5yds	24 x 23.5 x 28.5	57 x 33 x 26cm
	7.5cm x 5yds	24 x 23.5 x 41	57 x 33 x 26cm
	10.5cm x 5yds		57 x 33 x 26cm
	1.25cm x 10m		57 x 33 x 26cm
	2.5cm x 10m		
	5cm x 10m		
	7.5cm x 10m		
	10cm x 10m		



NON - WOVEN

CAP

- Material: PP non woven/SMS
- Weight: 20 gsm, 25gsm, 30gsm etc.
- Type: with tie or elastic
- Size: 62 x 12.5cm / 63 x 13.5cm
- Weight: Blue, green, yellow etc.
- Packing: 10pcs/bag, 100 pcs/ctn



COVERALL

- Material: PP/SMS/SF/MP
- Size: S,M,L,XL,XXL
- Weight: 35gsm, 40gsm, 50gsm,60gsm etc.
- With or whitout hood/shoe cover
- Elastic cuff or kitted cuff
- Color: White, blue, yellow etc.
- Packing: 1pc/pouch, 50pcs/ctn (Sterile)
- 5pcs/bag, 100 pcs/ctn (Non sterile)



AAMI SURGICAL GOWN

- Material: AAMI SMS
- Weight: 50g, 60g etc
- Size: S,M,L,XL,XXL
- Elastic cuff and knitted cuff
- Color: Blue, Light blue, etc.
- Packing: 1pc/blister, 25pcs/ctn (sterile)



TRANSVERSE INCISION

COMPONENTS	DESCRIPTION	SIZE	QUANTITY
Back Table cover	50g PE film	105 x 200cm	1pcs
Medium Drape	35g SMS	100 x 150 cm	2pcs
Mayo stand cover	PE film + 35g SMS	58 x 137cm	1pcs
Hand towel	66gsm white paper	40 x 40cm	5pcs
Drape utility with tape	35 g SMS+1cm	50 x 60cm	4pcs
Laparotomy Drape	35g SMS	200 x 230cm	1pcs
Reinforced Surgical Gown	43 gsm	127 x 150cm	2pcs

With adhesive tape

Antistaic SMMS



GENERAL PACK

COMPONENTS	DESCRIPTION	SIZE	QUANTITY
Mayo stand cover	42 gsm blue PE with 30gsm hydrophilic SPP non woven in middle, SPP width 60cm	80 x 145cm	1pc
Op-tape	68 gsm Souulance non woven with adhesive	10 x 50cm	2pcs
Hand Towel	66gsm white paper	40 x 40cm	2pcs
Side drape	27 gsm PE fillm with 30gsm hydrophilic non woven 27gsm	75 x 90cm	2pcs
Head drape	27gsm PE film with 30gsm hydrophilic 60gSPP non woven	150 x 240cm	1pc
Foot drape	27gsm PE film with 30gsm hydrophilic SPP non woven	150 x 180cm	1pc
Reinforced Surgical Gown	43gsm Antistaitc SMMS	127 x 150cm	2pcs
Instrument Table Cover	42gsm PE film with 30gsm hydrophilic SPP non woven in middle, SPP width 60cm	150 x 200cm	1pc
Wrapper	35gsm Antistatic SMS	100 x 100cm	1pc



DIAPER

BABY DIAPER



ITEM #	L x W	SAP WEIGHT	WEIGHT/PC	PACKING	CARTON
	380 X 280	105 x 200cm	4g	52pcs/bag, 4bags/ctn	46x25x30cm
	450 X 320	100 x 150 cm	5g	52pcs/bag, 4bags/ctn	46x25x30cm
	500 X 320	58 x 137cm	6g	52pcs/bag, 4bags/ctn	46x25x30cm
	540 X 320	40 x 40cm	7g	52pcs/bag, 4bags/ctn	46x25x30cm



ADULT DIAPER

Standard type: anti-leakage PE film, leg elastics, Left/right tapes, frontal tape, leg cuffs

ITEM #	L x W	SAP WEIGHT	WEIGHT/PC	PACKING	CARTON
	800 x 650	7.5g	85g	10pcs/bag, 10bags/ctn	86x24.5x30cm
	900 x 750	9g	95g	10pcs/bag, 10bags/ctn	86x27.5x40cm
	980 x 800	10g	105g	10pcs/bag, 10bags/ctn	86x28.5x41cm

Standard type: anti-leakage PE film, leg elastics, Left/right tapes, frontal tape, leg cuffs, wetness indicator

ITEM #	L x W	SAP WEIGHT	WEIGHT/PC	PACKING	CARTON
	800 x 650	7.5g	85g	10pcs/bag, 10bags/ctn	86x24.5x30cm
	900 x 750	9g	95g	10pcs/bag, 10bags/ctn	86x27.5x40cm
	980 x 800	10g	105g	10pcs/bag, 10bags/ctn	86x28.5x41cm

INFUSION SET

- Medical grade PVC, EO sterile
- Leuer lock or luer slip
- With or without needle
- Drip chamber with or without filter
- Spike with or without air vent
- Y Type, flash bulb injection, ordinary type, burette type
- Tube length: 1.35m, 1.5m, 2m or customized
- PE or blister packing

IV CANNULA

- Medical grade PVC + stainless steel
- Pen type, with wings, with injection port, Y type
- 14G, 16G, 18G, 20G, 22G, 24G, 26G
- Orange, grey, green, pink, blue, yellow, violet
- Sterilized by EO gas, non-toxic, non-pyrogenic
- Blister pack

COMPONENTS	DESCRIPTION	SIZE	QUANTITY
Back Table cover	50g PE film	105 x 200cm	1pcs
Medium Drape	35g SMS	100 x 150 cm	2pcs
Mayo stand cover	PE film + 35g SMS	58 x 137cm	1pcs
Hand towel	66gsm white paper	40 x 40cm	5pcs
Drape utility with tape	35 g SMS+1cm	50 x 60cm	4pcs
Laparotomy Drape	35g SMS	200 x 230cm	1pcs
Reinforced Surgical Gown	43 gsm	127 x 150cm	2pcs



NEEDLE FOR ANESTHESIA

- Medical grade PVC + stainless steel
- Type: Quincke tip or pencil-point
- 16G, 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
- Length: 45mm, 60mm, 75mm, 90mm, 120mm, 150mm etc
- With or without introducer available
- Sterilized by EO gas, sterile, non-toxic, non-pyrogenic




ITEM #	SIZE	INTRODUCER	COLOR
	18G x 90mm	None	Pink
	19G x 90mm	None	Cream
	20G x 90mm	None	Yellow
	21G x 90mm	None	Green
	22G x 90mm	None	Black
	23G x 90mm	19G x 38mm	Blue
	24G x 90mm	20G x 38mm	Purple
	25G x 90mm	20G x 38mm	Orange
	26G x 90mm	20 G x 38mm	Brown
	27G x 90mm	22G x 38mm	Grey

ITEM #	SIZE	INTRODUCER	COLOR
	22G x 90mm	None	Black
	23G x 90mm	19G x 38mm	Blue
	24G x 90mm	19G x 38mm	Purple
	25G x 90mm	20G x 38mm	Orange
	26G x 90mm	20G x 38mm	Brown
	27G x 90mm	22 G x 38mm	Medium Gray

FEEDING TUBE

- Medical grade PVC
- Transparent or frosted
- Tube length: 40cm + 2cm (Infant), 120 cm + 2cm (Adult)
- Size: F4, F5, F6, F8, F10, F12, F14, F16, F18, F20, F22.
- PE Packing or blister packing
- Sterilized by EO GAS



NEBULIZER MASK

- Medical grade PVC
- Adjustable nasal clip
- Size: S/ M / L/XL
- Suitable for adult, pediatric and infant
- Color: Green or white
- Tube length: 2m, 6m, OD



BREATHING CIRCUIT

- Material: Medical grade PVC
- Normal length: 2.0m, 2.4m, can be customized
- Model: Adult, child
- Color: Transparent, green, white
- Diameter: 22mm, 15mm
- Length: 1.6m / 1.8m / 2m / 3m
- Application: Anesthesia operation



ITEM #	SPECIFICATION	PACKING	CARTON SIZE
	22mm	100 pcs / ctn	74 x 39 x 35 (5cm)
	15mm	100 pcs / ctn	74 x 39 x 35 (5cm)

BLOOD BAG

- Main material: Medical grade PVC
- Bag size choice: 200, 250, 350, 400, 450, 500ml
- Type: Single, double, triple
- Material Origin: China or Netherland
- Function: Blood container
- Color: Transparent



ITEM #	ANTI-COAGULANT	PACKING	CARTON SIZE
	22mm	100 pcs / ctn	74 x 39 x 35 (5cm)
	15mm	100 pcs / ctn	74 x 39 x 35 (5cm)

ADULT URINE BAG

- Material: Made of PVC, non-toxic
- Size: 500ml, 750ml, 1000ml, 1500ml, 2000ml
- Thickness: 1.2mm
- Medical PVC tube: 0.65 x 90cm
- Milky or transparent
- Supplied in individual peelable polybag or blister pack sterile
- Type: Bottom with T valve / screw valve / without valve (Single tube)
- Sterilized by EO gas

ITEM #	SPECIFICATION	PACKING	SPECIFICATION
	2000ml w/o bottom outlet	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 29cm
	2000 ml with pull-push valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 29cm
	2000ml with screw valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 31cm
	200ml with T valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 33cm



AUTOCLAVE STERILIZATION INDICATOR TAPE

- Innovate designs are per request
- Sterilization method: STEAM and ETO
- Sizes: 19mm x 50m
- Material: Masking tape

STERILIZATION INDICATOR TAPE

ITEM #	SPECIFICATION	PACKING	SPECIFICATION
	2000ml w/o bottom outlet	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 29cm
	2000 ml with pull-push valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 29cm
	2000ml with screw valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 31cm
	200ml with T valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 33cm



ITEM #	SPECIFICATION	PACKING	SPECIFICATION
	2000ml w/o bottom outlet	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 29cm
	2000 ml with pull-push valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 29cm
	2000ml with screw valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 31cm
	200ml with T valve	1pc/ individual pack, 10pcs/PE bag, 25bags/ctn	55 x 41 x 33cm



HNM MEDICAL



20855 NE 16 Ave
Suite C15, Miami,
FL 33179



1.866.291.8498
1.888.608.9788



www.hnmmedical.com
info@hnmmedical.com