


HNM Tag-it is our new line of information tags developed to help medical professionals identify and document a wide range of instructions. Whether you need to identify instruments, trays or carts, both before or after sterilization, HNM Tag-It has got you covered, these durable and tear-resistant tags are perfect for indication special handling instructions.

WHY YOU'LL LOVE THEM:

- * The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2 F (134 C).
- * The synthetic paper is durable and will not wrinkle or absorb water. Made from recyclable plastic promoting environmental sustainability.
- * Robust and easy to use tie wraps designed to securely attach to a wide range of trays, instruments, and other medical devices.
- * Available in over 8 languages.



All tags available with indicator dots
Provides visual confirmation that the tag was processed



KEY PRODUCT FEATURES

- * High Quality Permanent Marker
- * Steam or Plasm Indicator Dots provide confirmation of processing.
- * Eco-friendly, Waterproof & Tearproof Tags
- * Large Quantities: Each bundle of HNM Tag-It comes with 100 versatile tags, to quickly identify instruments that require attention.
- * Available in multiple linguistic choices.

| ITEM# | COLOR | DESCRIPTION | SIZE | PACKAGING |
|--------------|------------|------------------------|--------------------------------|-----------|
| HNM-TI 3502 | Blue | VENDOR | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 3503 | Green | LOANER | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 3504 | Red | PRIORITY | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 3505 | Yellow | RE-REUN | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 3506 | Pink | MISSING | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 3507 | Orange | BROKEN | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 3509 | Black | BLANK | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 35B1 | White | BLANK | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 35B4 | Red | BROKEN | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 35C22 | Light Blue | VENDOR CHECK IN + OUT | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 35L3 | Green | LOANER CHECK IN + OUT | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI 35T2 | Blue | QUICK TURNOVER | 3" X 5" / 8 cm X 13 cm | 100 / PKG |
| HNM-TI PM09 | Black | FINE PERMANENT MARKERS | 1.5" X 5.5" / 3.81 cm X 12.7cm | 36 / Bag |

All information tags are single use. PACKING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker

PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for stearn dots.

PROCESSING INDICATOR DOTS: Add suffix "PL" to end of code for stearn dots.

infection control & sterile processing



VENDOR / BLUE

ITEM # **SIZE**
 HNM-TI 3502 3" X 5"
 (8 cm X 12 cm)

VENDOR

Company Name _____

Tray Name _____

Date of Surgery _____

TO REORDER TAG-IT™: 1.800.661.5432 WWW.BATRIK.COM TI-3502Z

LOANER / GREEN

ITEM # **SIZE**
 HNM-TI 3503 3" X 5"
 (8 cm X 12 cm)

LOANER

Loaned From _____

Loaned To _____

By _____ Date _____

TO REORDER TAG-IT™: 1.800.661.5432 WWW.BATRIK.COM TI-3503

PRIORITY / RED

ITEM # **SIZE**
 HNM-TI 3504 3" X 5"
 (8 cm X 12 cm)

PRIORITY

Room # _____

Time Needed _____

Initials _____

TO REORDER TAG-IT™: 1.800.661.5432 WWW.BATRIK.COM TI-3504

RE-RUN / YELLOW

ITEM # **SIZE**
 HNM-TI 3505 3" X 5"
 (8 cm X 12 cm)

RE-RUN

Room # _____

TO REORDER TAG-IT™: 1.800.661.5432 WWW.BATRIK.COM TI-3505

MISSING / PINK

ITEM # **SIZE**
 HNM-TI 3506 3" X 5"
 (8 cm X 12 cm)

MISSING

TO REORDER TAG-IT™: 1.800.661.5432 WWW.BATRIK.COM TI-3506

BROKEN / ORANGE

ITEM # **SIZE**
 HNM-TI 3507 3" X 5"
 (8 cm X 12 cm)

BROKEN

DO NOT USE

DATE: _____ PHONE: _____

DEPARTMENT/CONTACT: _____

ITEM: _____

CRITICAL: (☒) Yes No

TO REORDER TAG-IT™: 1.800.661.5432 WWW.BATRIK.COM TI-3507

BLANK / WHITE

ITEM # **SIZE**
 HNM-TI 35B1 3" X 5"
 (8 cm X 12 cm)

TO REORDER TAG-IT™: 1.800.661.5432 BATRIK MEDICAL MANUFACTURING INC. TI-35B1

BROKEN / RED

ITEM # **SIZE**
 HNM-TI 35B4 3" X 5"
 (8 cm X 12 cm)

BROKEN DO NOT USE

DATE: _____ PHONE: _____

DEPARTMENT/CONTACT: _____

ITEM: _____

CRITICAL: (☒) Yes No

INSTRUMENT SET: _____

ASSEMBLED BY: _____

TI-35B4

VENDOR IN + OUT / LIGHT BLUE

ITEM # **SIZE**
 HNM-TI 35C22 3" X 5"
 (8 cm X 12 cm)

VENDOR

Tray Received From: _____

Tray Name: _____

Received By: _____

Date: _____ Surgery Date: _____

CPD Notified Prior to Receipt of Trays
 Received >24 Hours Before Case
 Inventory List Available
 Inventory Check Complete
 All Instruments in Good Condition
 Containers in Good Condition
 All Trays Labeled (Patient & Surgeon)
 All Trays <25 pounds
 Written Instructions For Use Available
 Disassembled/Decontaminated

WWW.BATRIK.COM TI-35C22

VENDOR

Check Out

Disassembled/Decontaminated
 Inventory Check Complete
 All Items Complete and Accounted For

Returned By: _____

Received By: _____

Date: _____

Notes: _____

TO REORDER TAG-IT™: 1.800.661.5432 TI-35C22

LOANER IN + OUT / GREEN

LOANER

Check In

Tray Received From: _____

Tray Name: _____

Received By: _____

Date: _____ Surgery Date: _____

CPD Notified Prior to Receipt of Trays
 Received >24 Hours Before Case
 Inventory List Available
 Inventory Check Complete
 All Instruments in Good Condition
 Containers in Good Condition
 All Trays Labeled (Patient & Surgeon)
 All Trays <25 pounds
 Written Instructions For Use Available
 Disassembled/Decontaminated

WWW.BATRIK.COM TI-35L3

LOANER

Check Out

Disassembled/Decontaminated
 Inventory Check Complete
 All Items Complete and Accounted For

Returned By: _____

Received By: _____

Date: _____

Notes: _____

TO REORDER TAG-IT™: 1.800.661.5432 TI-35L3

ITEM # **SIZE**
 HNM-TI 35L3 3" X 5"
 (8 cm X 12 cm)

QUICK TURNOVER / BLUE

QUICK TURNOVER

Tech Initials: _____

Room Number: _____

Item Number: _____

Time Needed: _____

Decontam/LLD/ILD/HLDD
 Immediate Use Sterilization
 Terminal Sterilization in SPD

TO REORDER TAG-IT™: 1.800.661.5432 WWW.BATRIK.COM TI-35T2

ITEM # **SIZE**
 HNM-TI 35T2 3" X 5"
 (8 cm X 12 cm)