

HNM Scope Tag-Them are tear resistant, double sided tags that will prove indispensable throughout the reprocessing journey.

Our tags make it easy to record important information, including processing dates, with a user friendly feature that allows you to select the month, day, and date.

Crafted from waterproof synthetic paper, these tags feature a self looping design, making it effortless to attach them to the scope without the need for additional fastening materials.



KEY FEATURES

- The paper, tie wraps, and smudge proof ink are designed to withstand high temperatures of up to 273.2° F (134°C) and are compatible with Gas plasma sterilization STERRAD®100 NXSterilization System).
- Tags include the month, day, and date selection on one side, making it convenient to mark the processing date accurately.

COLOR-CODED AND CUSTOMIZABLE

- Available in seven colors
- The tags facilitate color coding for easy identification.

PRODUCT SELECTION

ITEM #	COLOR	SIZE	PACKAGING
HNM-ST05	Yellow	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
HNM-ST02	Blue	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
HNM-ST03	Green	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
HNM-ST07	Orange	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
HNM-ST06	Purple	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
HNM-ST04	Red	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
HNM-ST01	White	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG

All information tags are single use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.

SELF LOOPING

- The self looping design allows for easy attachment to the scope.
- Eliminating the need for additional fastening materials.

SINGLE-USE AND LATEX FREE

- HNM Scope Tag-Them are designed for single use and are latex free.
- Promoting hygiene and reducing the risk of cross contamination.

