



# Single-use **Laparoscopic** **Aspirating Needle**



- ✓ Stainless steel cannula with 18g needle for laparoscopic aspiration
- ✓ Extra long 335mm length suitable for all patients of any BMI
- ✓ Available in 3.5 and 5mm diameter cannulas complete with Luer connection

**1**

No burden of sterilization costs  
No repair costs for bent or damaged needles  
No shortage of instruments while undergoing repair

**2**

A delicate Single-use instrument that is sharp every time

**3**

Single-use eliminates the risk of cross contamination

